

Anti-CD49C (Integrin alpha3)

Mouse IgG1, kappa
P1B5

BioLegend
From Revvity



100 µL

Artikel-Nr.: B921601 | **BioLegend** | Herstellernummer: 921601

153,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: FC, IHC-F, IP
Clone: P1B5
Format: ASCITES
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 µl

Anti-CD49C (Integrin alpha3)

Clone: P1B5
Isotype: Mouse IgG1, kappa
Other_Names: GAP-B3, GAPB3, ILNEB, MSK18, VCA-2, VL3A, VLA3a, FRP-2, galactoprotein B3, VLA-3 subunit alpha, CD49 antigen-like family member C, antigen identified by monoclonal antibody J143, very late activation protein 3 receptor, alpha-3 subunit
Size: 100 µL

Reactivity: Human
Format: ASCITES
App_Abbrev: IHC-F, FC, IP
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11033>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>