

GoInVivo Purified anti-human CD366 (Tim-3)

Mouse IgG1, kappa
F38-2E2

BioLegend
From Revvity



500 mg

Artikel-Nr.: B345036 | BioLegend | Herstellernummer: 345036

18.964,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Block, Costim, FC

Clone: F38-2E2

Format: GOINVIVO

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 500 mg

GoInVivo Purified anti-human CD366 (Tim-3)

Clone: F38-2E2

Isotype: Mouse IgG1, kappa

Other_Names: T cell immunoglobulin and mucin domain containing protein 3, hepatitis virus cellular receptor 2, CD366

Size: 500 mg

Reactivity: Human

Format: GOINVIVO

App_Abbrev: FC, Costim, Block

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12235>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>