

Purified anti-Arginase 1

Mouse IgG1, kappa
O94E6/ARG1

BioLegend
From Revvity



100 µg

Artikel-Nr.: B678802 | **BioLegend** | Herstellernummer: 678802

203,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: IHC-P, KO/KD-WB, WB

Clone: O94E6/ARG1

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-Arginase 1

Clone: O94E6/ARG1

Isotype: Mouse IgG1, kappa

Other_Names: Arginase 1, liver type Arginase, Type I Arginase, ARG1

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: WB, KO/KD-WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12106>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>