

Purified anti-human LIF

Goat Polyclonal IgG

Poly5233

BioLegend
From Revvity



100 µg

Artikel-Nr.: B523301 | **BioLegend** | Herstellernummer: 523301

210,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ELISA Capture, WB

Clone: Poly5233

Format: PURE

Isotype: Goat Polyclonal IgG

Species of Origin: Goat

Verpackung: 100 µg

Purified anti-human LIF

Clone: Poly5233

Isotype: Goat Polyclonal IgG

Other_Names: Leukemia Inhibitory Factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture, WB

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12804>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>