

Purified anti-human Adiponectin (Acrp30)

Mouse IgG2a, kappa
62-9

BioLegend
From Revvity



500 µg

Artikel-Nr.: B528702 | [BioLegend](#) | Herstellernummer: 528702

306,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Capture

Clone: 62-9

Format: PURE

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 500 µg

Purified anti-human Adiponectin (Acrp30)

Clone: 62-9

Isotype: Mouse IgG2a, kappa

Other_Names: GBP-28, apM1, AdipoQ, and Acrp30

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12729>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>