

Purified anti-human Arginase I

Mouse IgG2b, kappa
14D2C43

BioLegend
From Revvity



100 µg

Artikel-Nr.: B369702 | **BioLegend** | Herstellernummer: 369702

203,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC

Clone: 14D2C43

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human Arginase I

Clone: 14D2C43

Isotype: Mouse IgG2b, kappa

Other_Names: ARG1, Liver Arginase, Type 1 Arginase, Arginase-1, Liver-Type Arginase

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12911>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>