

Purified anti-complement C3b/iC3b

Mouse IgG1, kappa
3E7/C3b

BioLegend
From Revvity



100 µg

Artikel-Nr.: B846102 | **BioLegend** | Herstellernummer: 846102

138,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Block, Direct ELISA, ELISA Detection, FC, ICC, RIA

Clone: 3E7/C3b

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-complement C3b/iC3b

Clone: 3E7/C3b

Isotype: Mouse IgG1, kappa

Other_Names: Complement Component 3, C3b, iC3b, C3bi

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, Direct ELISA, Block, ELISA Detection, ICC, RIA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12958>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>