

Purified anti-complement C3/C3b/iC3b

Mouse IgG1, kappa
6C9/C3b

BioLegend
From Revvity



100 µg

Artikel-Nr.: B846502 | **BioLegend** | Herstellernummer: 846502

159,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA, ELISA Detection, FC, WB

Clone: 6C9/C3b

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-complement C3/C3b/iC3b

Clone: 6C9/C3b

Isotype: Mouse IgG1, kappa

Other_Names: Complement Component 3, C3b, iC3b, C3bi

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, FC, Direct ELISA, ELISA Detection

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=12536>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>