

Purified anti-RNA Polymerase II RPB1

Mouse IgM, kappa
H14

BioLegend
From Revvity



100 µg

Artikel-Nr.: B920304 | **BioLegend** | Herstellernummer: 920304

194,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ICC, IP, WB
Clone: H14
Format: PURE
Isotype: kappa, Mouse IgM
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-RNA Polymerase II RPB1

Clone: H14
Isotype: Mouse IgM, kappa
Other_Names: DNA-directed RNA polymerase II subunit RPB1, RNA polymerase II subunit B1, DNA-directed RNA polymerase II subunit A, RNA-directed RNA polymerase II subunit RPB1, DNA-directed RNA polymerase III largest subunit
Size: 100 µg

Reactivity: Human, Mouse, Rat
Format: PURE
App_Abbrev: WB, ICC, IP
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=13076>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>