

APC anti-mouse CD96 (TACTILE)

Rat IgG1, kappa
3.3

BioLegend
From Revvity



25 µg

Artikel-Nr.: B131711 | [BioLegend](#) | Herstellernummer: 131711

120,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: 3.3
Format: APC
Isotype: kappa, Rat IgG1
Species of Origin: Rat
Verpackung: 25 µg

APC anti-mouse CD96 (TACTILE)

Clone: 3.3
Isotype: Rat IgG1, kappa
Other_Names: Tactile (T cell-activated increased late expression)
Size: 25 µg

Reactivity: Mouse
Format: APC
App_Abbrev: FC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=13998>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>