

MAGic Clamp COMBO

for 125 ml, 250 ml, 500 ml flasks



55H1000-MR-CMB Adjustable flask clamp, accepts 125ml, 250ml and 500ml Erlenmeyer flasks. (For a video demonstration, scan the QR code above.)

1 Stück

Artikel-Nr.: 55H1000-MR-CMB | **Benchmark** | Herstellernummer: H1000-MR-CMB

58,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Magic Clamp: ja

Verpackung: 1 Stück