

Biotin anti-human CD326 (EpCAM)

Mouse IgG2a, kappa
CO17-1A

BioLegend
From Revvity



100 µg

Artikel-Nr.: B369812 | **BioLegend** | Herstellernummer: 369812

212,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: CO17-1A
Format: BIOTIN
Isotype: kappa, Mouse IgG2a
Species of Origin: Mouse
Verpackung: 100 µg

Biotin anti-human CD326 (EpCAM)

Clone: CO17-1A
Isotype: Mouse IgG2a, kappa
Other_Names: Ep-CAM, tumor associated calcium signal transducer 1, epithelial cell surface antigen, epithelial glycoprotein 2, EGP2, adenocarcinoma associated antigen, TROP1
Size: 100 µg

Reactivity: Human
Format: BIOTIN
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14213>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>