

APC anti-human/mouse Granzyme B Recombinant

Mouse IgG1, kappa
QA16A02

BioLegend
From Revvity



100 tests

Artikel-Nr.: B372204 | **BioLegend** | Herstellernummer: 372204

261,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC

Clone: QA16A02

Format: APC

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 Tests

APC anti-human/mouse Granzyme B Recombinant

Clone: QA16A02

Isotype: Mouse IgG1, kappa

Other_Names: Granzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1, GZMB, CCP1, Asp-aseGranzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1, GZMB, CCP1, Asp-ase

Size: 100 tests

Reactivity: Human, Mouse

Format: APC

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14429>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>