

## Purified anti-human CD312 (EMR2)

Rat IgG2a, kappa  
W15101A

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B373102 | BioLegend | Herstellernummer: 373102

**244,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FC, IHC-F

**Clone:** W15101A

**Format:** PURE

**Isotype:** kappa, Rat IgG2a

**Species of Origin:** Rat

**Verpackung:** 100 µg

Purified anti-human CD312 (EMR2)

Clone: W15101A

Isotype: Rat IgG2a, kappa

Other\_Names: Epidermal growth factor module-containing mucin-like receptor 2, EMR2, AGRE2, ADGRE2

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, IHC-F

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14147>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>