

Purified anti-Nos2 (iNOS)

Rat IgG2b, kappa
W16030C

BioLegend
From Revvity



100 µg

Artikel-Nr.: B696802 | BioLegend | Herstellernummer: 696802

165,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ICC, ICFC, WB
Clone: W16030C
Format: PURE
Isotype: kappa, Rat IgG2b
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-Nos2 (iNOS)

Clone: W16030C
Isotype: Rat IgG2b, kappa
Other_Names: Nitric oxide synthase 2, nitric oxide synthase 2A, NOS2A, HEP-NOS, inducible nitric oxide synthase (iNOS), Nitric oxide synthase, inducible, Inducible NOS
Size: 100 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: WB, ICC, ICFC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14396>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>