

Recombinant Mouse S100A9 (carrier-free)

BioLegend
From Revvity



100 µg

Artikel-Nr.: B765406 | **BioLegend** | Herstellernummer: 765406

587,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: BA
Format: RECOM
Verpackung: 100 µg

Recombinant Mouse S100A9 (carrier-free)

Clone:

Isotype:

Other_Names: S100 Calcium-binding protein A9, Cystic Fibrosis Antigen B, Calgranulin B (CAGB, CGLB), Myeloid-related protein 14 (MRP14),

Size: 100 µg

Reactivity: Mouse

Format: RECOM

App_Abbrev: BA

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14560>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>