

## Purified anti-alpha-Synuclein, 80-96

Mouse IgG1, kappa  
A15115A

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B848302 | BioLegend | Herstellernummer: 848302

**175,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** Direct ELISA, IHC-P, WB

**Clone:** A15115A

**Format:** PURE

**Isotype:** kappa, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 100 µg

Purified anti-alpha-Synuclein, 80-96

Clone: A15115A

Isotype: Mouse IgG1, kappa

Other\_Names: Synuclein alpha-140, non-A4 component of amyloid, alpha-synuclein, isoform NACP140, non-amyloid beta component (NAP), PARK1, NACP, PD1

Size: 100 µg

Reactivity: Human, Mouse, Rat

Format: PURE

App\_Abbrev: IHC-P, Direct ELISA, WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=13977>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>