

Purified anti-alpha-Synuclein, 34-45

Mouse IgG1, kappa
A15110D

BioLegend
From Revvity



25 µg

Artikel-Nr.: B849101 | [BioLegend](#) | Herstellernummer: 849101

54,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: Direct ELISA, IHC-P
Clone: A15110D
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-alpha-Synuclein, 34-45

Clone: A15110D
Isotype: Mouse IgG1, kappa
Other_Names: Synuclein alpha-140, non-A4 component of amyloid, alpha-synuclein, isoform NACP140, non-amyloid beta component (NAP)
Size: 25 µg

Reactivity: Human
Format: PURE
App_Abbrev: IHC-P, Direct ELISA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14067>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>