

## Purified anti-alpha-Synuclein, 34-45

Mouse IgG1, kappa  
A15110D

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B849102 | **BioLegend** | Herstellernummer: 849102

**131,25 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Aktionsangebote:** Aktion  
**App\_Abbrev:** Direct ELISA, IHC-P  
**Clone:** A15110D  
**Format:** PURE  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-alpha-Synuclein, 34-45

Clone: A15110D  
Isotype: Mouse IgG1, kappa  
Other\_Names: Synuclein alpha-140, non-A4 component of amyloid, alpha-synuclein, isoform NACP140, non-amyloid beta component (NAP)  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: IHC-P, Direct ELISA  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14067>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>