

PE anti-human CD307c/FcRL3

Mouse IgG2b, kappa
H5/FcRL3

BioLegend
From Revvity



100 tests

Artikel-Nr.: B374406 | **BioLegend** | Herstellernummer: 374406

288,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: H5/FcRL3
Format: PE
Isotype: kappa, Mouse IgG2b
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD307c/FcRL3

Clone: H5/FcRL3
Isotype: Mouse IgG2b, kappa
Other_Names: Fc receptor-like 3 (FCRL3), Immune receptor translocation-associated protein 3 (IRTA3), SH2 domain-containing phosphatase anchor protein 2 (SPAP2), IFGP family protein 3 (IFGP3)
Size: 100 tests

Reactivity: Human
Format: PE
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14693>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>