

## HRP anti-alpha-Synuclein, 103-108

Mouse IgG1, kappa  
4B12/Synuclein

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B807806 | **BioLegend** | Herstellernummer: 807806

**190,50 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Aktionsangebote:** Aktion  
**App\_Abbrev:** IHC-P, WB  
**Clone:** 4B12/Synuclein  
**Format:** HRP  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

HRP anti-alpha-Synuclein, 103-108

Clone: 4B12/Synuclein  
Isotype: Mouse IgG1, kappa  
Other\_Names: NACP, PARK1, PARK4, PD1, Synuclein alpha-140, non-A4 component of amyloid, alpha-synuclein, isoform NACP140, non-A beta component of AD amyloid Parkinson disease  
Size: 100 µg

Reactivity: Human  
Format: HRP  
App\_Abbrev: IHC-P, WB  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14738>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>