


Goat anti-rabbit snd. antibody AzureSpectra 700

AzureSpectra 700 conjugated secondary antibody
for fluorescent Westerns



250 µl

Artikel-Nr.: 512128 | Azure

125,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Anregung: 690 nm

Anti-: Rabbit

Applikation: Fluoreszenz

Emission: 709 nm

Konjugat: Fluorophor

Origin: Goat

Verpackung: 250 µl