


Goat anti-mouse snd. antibody AzureSpectra 700

AzureSpectra 700 conjugated secondary antibody
for fluorescent Westerns



250 µl

Artikel-Nr.: 512129 | Azure

125,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Anregung: 690 nm
Anti-: Mouse
Applikation: Fluoreszenz
Emission: 709 nm
Konjugat: Fluorophor
Origin: Goat
Verpackung: 250 µl