

Purified anti-MRP-14 (S100A9)

Mouse IgG2b, kappa
A15105J

BioLegend
From Revvity



25 µg

Artikel-Nr.: B600301 | **BioLegend** | Herstellernummer: 600301

81,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: WB
Clone: A15105J
Format: PURE
Isotype: kappa, Mouse IgG2b
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-MRP-14 (S100A9)

Clone: A15105J
Isotype: Mouse IgG2b, kappa
Other_Names: S100A9, Migration inhibitory factor-related proteins-14, Calgranulin-B, Leukocyte L1 complex heavy chain, Calprotectin L1H subunit
Size: 25 µg

Reactivity: Human
Format: PURE
App_Abbrev: WB
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15363>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>