

Purified anti-BTK Phospho (Tyr551)/ITK Phospho...

Mouse IgG1, kappa
A16064A

BioLegend
From Revvity



25 µg

Artikel-Nr.: B646901 | [BioLegend](#) | Herstellernummer: 646901

90,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC, WB

Clone: A16064A

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-BTK Phospho (Tyr551)/ITK Phospho (Tyr512)

Clone: A16064A

Isotype: Mouse IgG1, kappa

Other_Names: BTK: Bruton's tyrosine kinase, Bruton agammaglobulinemia tyrosine kinase, AGMX1, ATK, IMD1, PSCTK1, XLA, BPK, ITK: Interleukin-2-inducible T-cell kinase, T-cell-specific kinase, EMT, LYK, PSCTK2

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: WB, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15335>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>