

Purified anti-Adenomatous polyposis coli (APC)

Mouse IgG1, kappa

Ali 12-28

BioLegend
From Revvity



25 µg

Artikel-Nr.: B850801 | **BioLegend** | Herstellernummer: 850801

84,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICC, WB

Clone: Ali 12-28

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-Adenomatous polyposis coli (APC)

Clone: Ali 12-28

Isotype: Mouse IgG1, kappa

Other_Names: APC, BTPS2, DP2, DP2.5, DP3, PPP1R46, WNT signaling pathway regulator, adenomatosis polyposis coli tumor suppressor, deleted in polyposis 2.5, protein phosphatase 1, regulatory subunit 46, truncated adenomatosis polyposis coli

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11214>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>