

Purified anti-Myelin Basic Protein

Mouse IgG1, kappa
P82H9

BioLegend
From Revvity



25 µg

Artikel-Nr.: B850901 | **BioLegend** | Herstellernummer: 850901

67,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: IHC-P, WB

Clone: P82H9

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-Myelin Basic Protein

Clone: P82H9

Isotype: Mouse IgG1, kappa

Other_Names: MBP, Golli-MBP, myelin A1 protein, myelin membrane encephalitogenic protein

Size: 25 µg

Reactivity: Human, Mouse, Rat

Format: PURE

App_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14946>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>