

Purified anti-DJ-1 (PARK7)

Mouse IgG2b, kappa
A16125E

BioLegend
From Revvity



25 µg

Artikel-Nr.: B851501 | **BioLegend** | Herstellernummer: 851501

85,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA, IHC-P, WB

Clone: A16125E

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-DJ-1 (PARK7)

Clone: A16125E

Isotype: Mouse IgG2b, kappa

Other_Names: Protein DJ-1, oncogene, DJ1, PARK7, Parkinson disease protein 7, Parkinsonism-Associated Deglycase, HEL-S-67p, epididymis secretory sperm binding protein Li 67p, GATD2, Protein Deglycase DJ-1

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, Direct ELISA, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15300>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>