

Purified anti-HMB-45

Mouse IgG1, kappa
HMB-45

BioLegend
From Revvity



100 µg

Artikel-Nr.: B911504 | **BioLegend** | Herstellernummer: 911504

249,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: IHC-F, IHC-P

Clone: HMB-45

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-HMB-45

Clone: HMB-45

Isotype: Mouse IgG1, kappa

Other_Names: Melanoma Marker antibody, Human Melanoma Black, Melanocyte protein PMEL, PMEL

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: IHC-P, IHC-F

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14906>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>