

Purified anti-Tata Binding Protein

Mouse IgG1, kappa
1TBP18

BioLegend
From Revvity



25 µg

Artikel-Nr.: B668305 | **BioLegend** | Herstellernummer: 668305

60,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: KO/KD-WB, WB
Clone: 1TBP18
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Tata Binding Protein

Clone: 1TBP18
Isotype: Mouse IgG1, kappa
Other_Names: TATA-box-binding protein, TATA sequence-binding protein, TATA-binding factor, TATA-box factor, Transcription initiation factor TFIID TBP subunit, TBP
Size: 25 µg

Reactivity: Human, Mouse
Format: PURE
App_Abbrev: WB, KO/KD-WB
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15151>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>