

## Ultra-LEAF Purified anti-CD95 (Fas)

Mouse IgG2b, kappa  
A16086F

**BioLegend**  
From Revvity



1 mg

Artikel-Nr.: B684402 | BioLegend | Herstellernummer: 684402

**572,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** Block

**Clone:** A16086F

**Format:** ULEAF

**Isotype:** kappa, Mouse IgG2b

**Species of Origin:** Mouse

**Verpackung:** 1 mg

Ultra-LEAF Purified anti-CD95 (Fas)

Clone: A16086F

Isotype: Mouse IgG2b, kappa

Other\_Names: Tumor necrosis factor receptor superfamily member 6 (TNFRSF6), ALPS1A, APO-1, Apoptosis antigen 1 (APT1), CD95, FAS1, FASTM, FS-7-associated surface antigen

Size: 1 mg

Reactivity: Human

Format: ULEAF

App\_Abbrev: Block

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=14941>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>