

Prime-Step Prestained Broad Range Protein Ladd...

BioLegend
From Revvity



50 µL

Artikel-Nr.: B773301 | **BioLegend** | Herstellernummer: 773301

21,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: WB

Format: ANC

Verpackung: 50 µl

Prime-Step Prestained Broad Range Protein Ladder

Clone:

Isotype:

Other_Names: protein ladder, protein marker, pre-stained protein ladder, pre-stained protein marker, protein molecular weight marker

Size: 50 µL

Reactivity:

Format: ANC

App_Abbrev: WB

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15579>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>