

AzureSpectra Labeling Kit 650

Kit for covalently labeling antibodies with a
fluorescent dye



1 Stück

Artikel-Nr.: 512187 | Azure

517,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Verpackung: 1 Stück

AzureSpectra Labeling Kit 650 - Kit for covalently labeling antibodies with a fluorescent dye – includes labeling buffer, dye solution, quenching solution, neutralization buffer and two G-25 columns – each kit includes sufficient reagents to perform two labelings of 50 micrograms of antibody