

## AzureSpectra Labeling Kit 700

Kit for covalently labeling antibodies with a  
fluorescent dye



1 Stück

Artikel-Nr.: 512188 | Azure

**517,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Verpackung:** 1 Stück

AzureSpectra Labeling Kit 700 - Kit for covalently labeling antibodies with a fluorescent dye – includes labeling buffer, dye solution, quenching solution, neutralization buffer and two G-25 columns – each kit includes sufficient reagents to perform two labelings of 50 micrograms of antibody