

Purified anti-Apo E, 154-162

Mouse IgG2b, kappa
A17065A

BioLegend
From Revvity



25 µg

Artikel-Nr.: B852801 | [BioLegend](#) | Herstellernummer: 852801

90,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA, IHC-P, WB

Clone: A17065A

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 25 µg

Purified anti-Apo E, 154-162

Clone: A17065A

Isotype: Mouse IgG2b, kappa

Other_Names: Apolipoprotein E, APO-E, Alzheimer Disease 2, LDLQC5, APOE, LPG, AD2

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, Direct ELISA, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=15693>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>