

Purified (azide-free) anti-Tau, 404-441

Mouse IgG1, kappa
Tau46

BioLegend
From Revvity



50 µL

Artikel-Nr.: B806604 | **BioLegend** | Herstellernummer: 806604

67,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: IHC-P, WB
Clone: Tau46
Format: PUREAF
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 50 µl

Purified (azide-free) anti-Tau, 404-441

Clone: Tau46

Isotype: Mouse IgG1, kappa

Other_Names: Microtubule-associated protein tau, PHF-tau, paired helical filament-tau, neurofibrillary tangle, microtubule-associated protein tau, isoform 4, G protein beta1/gamma2 subunit-interacting factor 1, DDPAC, FTDP-17, MAPTL, MSTD, MTBT1, MTBT2, PPND
Size: 50 µL

Reactivity: Human

Format: PUREAF

App_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11568>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>