

Biotin anti-mouse P2RY12

Rat IgG2b, kappa
S16007D

BioLegend
From Revvity



100 µg

Artikel-Nr.: B848007 | **BioLegend** | Herstellernummer: 848007

190,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: IHC-P

Clone: S16007D

Format: BIOTIN

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 100 µg

Biotin anti-mouse P2RY12

Clone: S16007D

Isotype: Rat IgG2b, kappa

Other_Names: P2Y12, ADPG-R, BDPLT8, HORK3, P2T(AC), P2Y(12)R, P2Y(AC), P2Y(ADP), P2Y(cyc), P2Y12P2Y12 platelet ADP receptor, purinergic receptor P2RY12, purinergic receptor P2Y, G-protein coupled, 12

Size: 100 µg

Reactivity: Mouse

Format: BIOTIN

App_Abbrev: IHC-P

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=16149>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>