

## Purified anti-Tollip

Mouse IgG1, kappa  
O92A4

**BioLegend**  
From Revvity



100 µg

Artikel-Nr.: B854402 | BioLegend | Herstellernummer: 854402

**221,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** WB  
**Clone:** O92A4  
**Format:** PURE  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-Tollip

Clone: O92A4  
Isotype: Mouse IgG1, kappa  
Other\_Names: Tollip, Toll Interacting Protein  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: WB  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=16407>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>