

Purified Non-muscle Myosin Heavy Chain II-A

Rabbit Polyclonal IgG

Poly19098

BioLegend
From Revvity



25 µL

Artikel-Nr.: B909802 | BioLegend | Herstellernummer: 909802

72,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, WB

Clone: Poly19098

Format: PURE

Isotype: Rabbit Polyclonal IgG

Species of Origin: Rabbit

Verpackung: 25 µL

Purified Non-muscle Myosin Heavy Chain II-A

Clone: Poly19098

Isotype: Rabbit Polyclonal IgG

Other_Names: Myosin-9, NMMHC-A, MYH9 variant protein, myosin heavy chain 9, non-muscle myosin heavy chain A, cellular myosin heavy chain, type A, non-muscle myosin heavy polypeptide 9

Size: 25 µL

Reactivity: Human, Mouse, Rat

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=11471>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>