

Purified anti-human CCL13 (MCP-4)

Mouse IgG1, kappa
3G4

BioLegend
From Revvity



500 µg

Artikel-Nr.: B537602 | **BioLegend** | Herstellernummer: 537602

348,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ELISA Capture
Clone: 3G4
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human CCL13 (MCP-4)

Clone: 3G4
Isotype: Mouse IgG1, kappa
Other_Names: Monocyte chemoattractant protein 4 (MCP4), CK?10, New CC chemokine 1 (NCC1), Small inducible cytokine subfamily a, member 13 (SCYA13)
Size: 500 µg

Reactivity: Human
Format: PURE
App_Abbrev: ELISA Capture
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17115>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>