

Purified anti-Arl13b

Mouse IgG2a, kappa
N295B/66

BioLegend
From Revvity



100 µg

Artikel-Nr.: B857602 | **BioLegend** | Herstellernummer: 857602

99,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: IHC-P, WB

Clone: N295B/66

Format: PURE

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-Arl13b

Clone: N295B/66

Isotype: Mouse IgG2a, kappa

Other_Names: ADP-ribosylation factor-like protein 13B, ADP-ribosylation factor-like protein 2-like 1, ARL2-like protein 1 and Arl2l1

Size: 100 µg

Reactivity: Mouse, Rat

Format: PURE

App_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=16812>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>