

Goat anti-rabbit second. antibody AzureSpectra 490

AzureSpectra 490 conjugated secondary antibody
for fluorescent Westerns



250 µl

Artikel-Nr.: 512206 | Azure

125,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Anregung: 491 nm

Anti-: Rabbit

Applikation: Fluoreszenz

Emission: 515 nm

Konjugat: Fluorophor

Origin: Goat

Verpackung: 250 µl