

TotalSeq-B0259 anti-human Hashtag 9

Mouse IgG1, kappa (all clones)

LNH-94; 2M2

BioLegend
From Revvity



10 µg

Artikel-Nr.: B394647 | **BioLegend** | Herstellernummer: 394647

310,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: PG

Clone: LNH-94; 2M2

Format: TOTALSEQ_B

Isotype: kappa (all clones), Mouse IgG1

Species of Origin: Mouse

Verpackung: 10 µg

TotalSeq-B0259 anti-human Hashtag 9

Clone: LNH-94; 2M2

Isotype: Mouse IgG1, kappa (all clones)

Other_Names: Na+/K+ ATPase beta3, ATP1B3, beta2M, beta2-M, beta2-microglobulin b2-M, b2M

Size: 10 µg

Reactivity: Human

Format: TOTALSEQ_B

App_Abbrev: PG

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17769>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>