

LEGEND MAX Human PTX3 ELISA Kit

BioLegend
From Revvity



1 Pre-coated Plate

Artikel-Nr.: B446807 | **BioLegend** | Herstellernummer: 446807

277,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ELISA

Format: LEGEND MAX

Species of Origin: N/A

Verpackung: 1 Pre-coated Plate

LEGEND MAX Human PTX3 ELISA Kit

Clone:

Isotype:

Other_Names: PTX-3, Pentraxin 3, TSG-14, TNFAIP5, Tumor Necrosis Factor Inducible Gene 14 Protein, Tumor Necrosis Factor Alpha-Induce Protein 5, TNF Alpha-Induced Protein 5

Size: 1 Pre-coated Plate

Reactivity: Human

Format: LEGEND MAX

App_Abbrev: ELISA

Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17758>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>