

Purified anti-SEN2

Rat IgG2b, kappa
W17126B

BioLegend
From Revvity



100 µg

Artikel-Nr.: B618402 | BioLegend | Herstellernummer: 618402

203,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: WB
Clone: W17126B
Format: PURE
Isotype: kappa, Rat IgG2b
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-SEN2

Clone: W17126B
Isotype: Rat IgG2b, kappa
Other_Names: SUMO Specific Peptidase 2, SUMO1/Sentrin/SMT3 Specific Peptidase 2, SUMO1/Sentrin/SMT3 Specific Protease 2, Sentrin/SUMO-Specific Protease SEN2, Sentrin (SUMO1) Specific Peptidase 2, SMT3-Specific Isopeptidase 2, Sentrin-Specific Protease 2
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: WB
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17737>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>