

Purified anti-Ataxin-2

Mouse IgG2b, kappa
A17203C

BioLegend
From Revvity



100 µg

Artikel-Nr.: B868002 | **BioLegend** | Herstellernummer: 868002

165,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: IHC-P, WB

Clone: A17203C

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-Ataxin-2

Clone: A17203C

Isotype: Mouse IgG2b, kappa

Other_Names: Spinocerebellar ataxia type 2 protein, Trinucleotide repeat-containing gene 13 protein, ATX2, SCA2, TNRC13, SCA2

Size: 100 µg

Reactivity: Human, Mouse, Rat

Format: PURE

App_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17947>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>