

Purified anti-Cortactin (CTTN)

Mouse IgG1, kappa
4F11/CTTN

BioLegend
From Revvity



25 µg

Artikel-Nr.: B868101 | [BioLegend](#) | Herstellernummer: 868101

85,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, WB
Clone: 4F11/CTTN
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Cortactin (CTTN)

Clone: 4F11/CTTN
Isotype: Mouse IgG1, kappa
Other_Names: Amplexin, Src Substrate Cortactin, Oncogene EMS1
Size: 25 µg

Reactivity: Human, Mouse, Rat
Format: PURE
App_Abbrev: WB, ICC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=17916>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>