

## Purified anti-Aldose Reductase (AKR1B1)

Mouse IgG1, kappa  
4G9-B6-F6

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B870301 | [BioLegend](#) | Herstellernummer: 870301

**85,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** IHC-P, WB

**Clone:** 4G9-B6-F6

**Format:** PURE

**Isotype:** kappa, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 25 µg

Purified anti-Aldose Reductase (AKR1B1)

Clone: 4G9-B6-F6

Isotype: Mouse IgG1, kappa

Other\_Names: Aldo-keto reductase family 1 member B1, Aldehyde reductase, Aldose reductase, AR

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App\_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=18678>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>