

Filterspitzen 20 µl BevelPoint

steril

Nachfolgeartikel für Art.-Nr. 701041



10 Racks à 96 Stück

Artikel-Nr.: LC2555 | Labcon | Herstellernummer: 1055-965-018-9

77,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Farbe: farblos

Filter: ja

Länge: normal

Low Binding: nein

Maximalvolumen: 20 µl

Pipettenspitzen Linie: Labcon

Sterilität: steril

Verpackung: 10 Racks à 96 Stück

Filterspitze passend für:

- Brand Transferpette 2-20 µl, 8/12 Kanal 2-20 µl, 2.5-25 µl
- Eppendorf Reference 2 2-20 µl
- Eppendorf Research Pro 2-20 µl
- Eppendorf Xplorer (plus) 2-20 µl
- Gilson Neo 2-20 µl
- Gilson Pipetman Classic G, L 2-20 µl
- Gilson Pipetman M 2-20 µl, 8/12 Kanal 1-20 µl
- Gilson Pipetman Ultra 2-20 µl, 8/12 Kanal 1-20 µl
- Nichiryo 20 µl
- Rainin Classic 2-20 µl
- Rainin EDP1 2-20 µl
- Rainin EDP3 Plus 2-20 µl
- Rainin Pipet Lite SL 2-20 µl
- Sartorius mLine 2-20 µl
- Sartorius Proline plus, Tacta 2-20 µl
- Socorex Acura 825, Manual XS 2-20 µl
- Thermo Digital, F1, F2 2-20 µl
- Thermo Focus 2-20 µl