

Purified anti-GIT1

Rat IgG2a, kappa
W18253B

BioLegend
From Revvity



25 µg

Artikel-Nr.: B935601 | **BioLegend** | Herstellernummer: 935601

81,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, IP, WB

Clone: W18253B

Format: PURE

Isotype: kappa, Rat IgG2a

Species of Origin: Rat

Verpackung: 25 µg

Purified anti-GIT1

Clone: W18253B

Isotype: Rat IgG2a, kappa

Other_Names: G Protein-Coupled Receptor Kinase Interacting, ArfGAP 1, GRK-Interacting Protein 1, ARF GAP GIT1, CAT1, Cool-associated and tyrosine phosphorylated protein 1, ARF GTPase-activating protein GIT1

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: WB, ICC, IP

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=19047>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>